



STEP BY STEP PROCEDURE TO DOWNLOAD THE LabVIEW RUN TIME ENGINE

1. Click on the link LABVIEW RUNTIME ENGINE



2. Click on LVRTE90std.exe to download the exe file



New user has to make an account to download the Run Time Engine. The existing users having account in NI then can directly login to download the software.

https://sine.ni.com - Drivers and Updates - National Instruments - Microsoft Internet Explorer

File Edit View Favorites Tools Help

NATIONAL INSTRUMENTS

My Profile NI Home

Drivers and Updates

Please complete the form below to continue or log into your existing profile on the right.

Benefits of Creating a Profile
Your free NI.com profile gives you unlimited access to ni.com features only available to registered users..

- Access MyNI - your personalized web page at ni.com
- Place and view online orders and quotes
- Download and view technical documents, evaluation

Login
E-Mail Address

Password

[Can't remember your password?](#)
[Continue >](#)

Make a new account here by filling the form OR Login if you have an account

* - indicates required field

Region

* Choose Your Country

* Language

My Information

* First Name
* Last Name
* E-Mail Address
* Password
* Confirm Password
* Area * Phone Number Ext.

Use this selection to tailor your ni.com experience to your region. As a worldwide company, NI offers localized product and application information.

Create your ni.com user ID by entering your name, e-mail address, and unique password. Use your e-mail address and password to login to your personalized ni.com account and receive the latest NI information that fits your needs. Remember that your password is case-sensitive.

3. After login you will again get the second window now click on LVRTE90std.exe

http://joule.ni.com - LabVIEW Run-Time Engine 2009 - Windows 2000/Vista x64/Vista x86/XP - (32- - Microsoft Internet Explorer

File Edit View Favorites Tools Help

NATIONAL INSTRUMENTS

Cart | Help Search

MyNI | Contact NI

Support

Product-Specific Support
Drivers and Updates
Product Reference
KnowledgeBase
Services Resource Center
Troubleshooting
Support Utilities
Discussion Forums

Request Support
from an engineer

Products & Services **Solutions** **Support** **NI Developer Zone** **Academic** **Events** **Company**

NI Home > Support > Drivers and Updates Select Your Country

LabVIEW Run-Time Engine 2009 - Windows 2000/Vista x64/Vista x86/XP - (32-bit Standard RTE)

[Print Document](#) 15 ratings | 3.80 out of 5

LabVIEW 2009 Run-Time Engine (Standard) (32-bit) for Windows

Available Downloads:
Standard Download: [LVRTE90std.exe](#) (160 MB) **Click here**
Download Language: Chinese (Simplified); English; French; German; Japanese; Korean
Product Line: LabVIEW
Version: 2009
Release date: 08-2009
Software type: Run-Time
Operating system: Windows 2000; Windows Vista x64; Windows Vista x86; Windows XP


Description

This document contains the LabVIEW 2009 Run-Time Engine (Standard) (32-bit). This Run-Time Engine must be installed on any Windows system where you plan to run executables that you build with the LabVIEW Application Builder. This version of the Run-Time Engine may be used only to run executables that were built with LabVIEW 2009 (32-bit).


Additionally, the Run-Time Engine allows your browser to display VIs that are embedded in Web pages.

If you were attempting to view a VI through a browser and you were directed to this site by your Web browser, the VI will display correctly after you install the corresponding Run-Time Engine: 32-bit for 32-bit browsers and 64-bit for 64-bit browsers. NOTE: Windows Vista x64 has the 32-bit version of Internet Explorer set to be the default browser.


4. After clicking you will get the window shown below. Click “Save” to save the file.




5. After saving the file downloading will start itself.



6. Double click on the icon of LVRTE90std



7. Unzip the file by clicking the option “Unzip”.




8. After doing unzip you will get a window shown below. Click “OK” to proceed.




9. Click “Next” to install Run time Engine




10. Click “Next” to select the primary installation directory.




11. Select “NI LabVIEW Run-Time Engine 2009” and click “next” .




12. Click “Next”




13. A window shown below will appear.




14. Accept the terms and click “Next”




15. Click “Next” to begin installation.



16. Installation begins as shown in the window



17. A window shown below will appear after the installation is complete. Click “Finish”.



18. It is recommended to restart the computer before running this software.

